



ptfesutures.com

Golnit[®] PTFE Suture

Periodontal, Bone Graft
and Implant Surgery

The ideal suture for
dental surgery



KEY ADVANTAGES

- Completely non-absorbable, 100% pure PTFE
- Monofilament – prevents bacterial wicking
- Biologically inert and chemically nonreactive
- Decreased inflammation and faster wound repair
- Soft – comfortable for patients
- Great handling characteristics and knot security












APPLICATIONS

- Dental bone grafting
- Soft tissue grafting
- Implant surgery
- Periodontal surgery
- Ridge augmentation
- Guided tissue regeneration (GTR)

WHY GOLNIT® PTFE SUTURES?

- The largest variety of PTFE sutures ever offered in the market
- Custom-made sutures for dental and aesthetic surgeries
- Needles made of highest grade, 302 series stainless steel wire
- Excellent needle penetration and smooth surface

We carry many more specifications for USP 6/0, 5/0, 4/0, 3/0, 2/0. Please contact us for additional sizes and lengths.

NEEDLE DESCRIPTION	NEEDLE LENGTH	AVAILABLE USP SIZES	PRODUCT CODES
Reverse Cutting 3/8 Circle ▾	24mm 	2/0	fRD202468 or mRD202468
		3/0	fRD302463 or mRD302463
Reverse Cutting 3/8 Circle ▾	21mm 	2/0	fRD202168 or mRD202168
		3/0	fRD302163 or mRD302163
		4/0	fRD402158 or mRD402158
Reverse Cutting 3/8 Circle ▾	20mm 	2/0	fRD202068 or mRD202068
		3/0	fRD302063 or mRD302063
		4/0	fRD402058 or mRD402058
Reverse Cutting 3/8 Circle ▾	19mm 	2/0	fRD201968 or mRD201968
		3/0	fRD301958 or mRD301958
		4/0	fRD401958 or mRD401958
		5/0	fRD501948 or mRD501948
Reverse Cutting 3/8 Circle ▾	18mm 	2/0	fRD201863 or mRD201863
		3/0	fRD301858 or mRD301858
		4/0	fRD401853 or mRD401853
Reverse Cutting 3/8 Circle ▾	16mm 	3/0	fRD301658 or mRD301658
		4/0	fRD401648 or mRD401648
		5/0	fRD501643 or mRD501643
		6/0	fRD601643 or mRD601643
Reverse Cutting 3/8 Circle ▾	13mm 	4/0	fRD401343 or mRD401343
		5/0	fRD501343 or mRD401343
		6/0	fRD601338 or mRD601338
Taper Point PERIO 1/2 Circle ☐	13mm 	4/0	fTH401343 or mTH401343
		5/0	fTH501338 or mTH501338
		6/0	fTH601333 or mTH601333
Reverse Cutting 3/8 Circle ▾	12mm 	4/0	fRD401243 or mRD401243
		5/0	fRD501243 or mRD501243
		6/0	fRD601238 or mRD601238
Reverse Cutting 3/8 Circle ▾	9mm 	5/0	fRD500933 or mRD500933
		6/0	fRD600933 or mRD600933
Taper Point 3/8 Circle ☐	9mm 	5/0	fTD500933 or mTD500933
		6/0	fTD600928 or mTD600928