

ePTFE Suture Ordering Information

THE PERFECT COMBINATION OF SOPHISTICATION AND STRENGTH



a surgical
suture
company

GOLNIT


ePTFE Suture

Golnit® ePTFE sutures are highly popular in cardiovascular surgeries. Thanks to its design, the ePTFE suture is exceptionally useful during procedures where the bleeding is a risk factor:

- Needle-hole bleeding reduction: The unique suture design, based on Golnit's proprietary technology, allows achieving a needle to suture diameter ratio of 1.
- Next Generation Golnit® sutures with Premium Sharp needles: Needles remain sharp pass after pass, limiting trauma to tissues.
- Superior physical and biological characteristics: ePTFE sutures are completely biocompatible and foster a benign tissue response. The sutures are also easy to handle, are very soft and flexible.




Golnit® ePTFE sutures are suitable for a wide range of procedures that include but are not limited to:

- Mitral valve repair and replacement
- Vascular reconstruction and repair
- Ventral and inguinal hernia repair

		Approx USP Size	7/0	6/0	5/0	4/0	3/0	2/0	0	
		Thread Diameter, mm	0,091-0,108	0,109-0,167	0,168-0,245	0,246-0,306	0,307-0,421	0,422-0,517	0,518-0,625	
		Gore Reference Thread Size:	CV-8	CV-7	CV-6	CV-5	CV-4	CV-3	CV-2	
		Needle Length, mm	Gore Reference Needle Code	Thread Length	Golnit Product Code (Gore Reference)					
Taper Point 3/8 Circle 	Double Armed 	6,4	TTc-7	7 cm	FTDVS806..d-7 (Gore Reference: 8D10)					
			TTc-7	61 cm	FTDVS806..d-61 (Gore Reference: 8K10)					
		9	TTc-9	7 cm	FTDVS809..d-7 (Gore Reference: 8D14)					
			TTc-9	46 cm	FTDVS809..d-46 (Gore Reference: 8J02)	FTDVS709..d-46 (Gore Reference: 7J02)	FTDVS609..d-46 (Gore Reference: 6J02)			
			TTc-9	61 cm	FTDVS809..d-61 (Gore Reference: 8K02)	FTDVS709..d-61 (Gore Reference: 7K02)	FTDVS609..d-61 (Gore Reference: 6K02)			
			TTc-9	76 cm	FTDVS809..d-76 (Gore Reference: 8M02)	FTDVS709..d-76 (Gore Reference: 7M02)	FTDVS609..d-76 (Gore Reference: 6M02)	FTDVS509..d-76		
			TTc-9	91 cm			FTDVS609..d-91			
		12	TTc-12	61 cm			FTDVS612..d-61 (Gore Reference: 6K10)			
			TTc-12	76 cm	FTDVS812..d-76 (Gore Reference: 8M04)	FTDVS712..d-76 (Gore Reference: 7M12)	FTDVS612..d-76 (Gore Reference: 6M10)			
		13	TTc-13	46 cm		FTDVS713..d-46 (Gore Reference: 7J06)	FTDVS613..d-46 (Gore Reference: 6J04)			
			TTc-13	61 cm		FTDVS713..d-61 (Gore Reference: 7K06)	FTDVS613..d-61 (Gore Reference: 6K04)	FTDVS513..d-61 (Gore Reference: 5K02)		
			TTc-13	76 cm	FTDVS813..d-76 (Gore Reference: 8M04)	FTDVS713..d-76 (Gore Reference: 7M06)	FTDVS613..d-76 (Gore Reference: 6M04)	FTDVS513..d-76 (Gore Reference: 5M02)	FTDVS413..d-76	
			TTc-13	91 cm			FTDVS613..d-91 (Gore Reference: 6N04)	FTDVS513..d-91 (Gore Reference: 5N02)	FTDVS413..d-91	
		17	TT-17	76 cm			FTDVS617..d-76	FTDVS517..d-76	FTDVS417..d-76	FTDVS317..d-76
			TT-17	91 cm				FTDVS517..d-91 (Gore Reference: 5N06)	FTDVS417..d-91 (Gore Reference: 4N04)	FTDVS317..d-91 (Gore Reference: 3N02)

Approx USP Size	7/0	6/0	5/0	4/0	3/0	2/0	0
Thread Diameter, mm	0,091-0,108	0,109-0,167	0,168-0,245	0,246-0,306	0,307-0,421	0,422-0,517	0,518-0,625
Gore Reference Thread Size:	CV-8	CV-7	CV-6	CV-5	CV-4	CV-3	CV-2

Taper Point 1/2 Circle	Double Armed	Needle Length, mm	Gore Reference Needle Code	Thread Length	Golnit Product Code (Gore Reference)							
		9	THc-9	76 cm	FTHVS809..d-76	FTHVS709..d-76	FTHVS609..d-76 (Gore Reference: 6M24)					
		13	THc-13	46 cm		FTHVS713..d-46 (Gore Reference: 7J10)	FTHVS613..d-46 (Gore Reference: 6J12)					
			THc-13	61 cm			FTHVS613..d-61 (Gore Reference: 6K12)	FTHVS513..d-61 (Gore Reference: 5K12)				
			THc-13	76 cm		FTHVS713..d-76 (Gore Reference: 7M10)	FTHVS613..d-76 (Gore Reference: 6M12)	FTHVS513..d-76 (Gore Reference: 5M12)	FTHVS413..d-76			
			THc-13	91 cm			FTHVS613..d-91	FTHVS513..d-91 (Gore Reference: 5N12)				
		18	TH-18	122 cm					FTHVS418..-122 (Gore Reference: 4U02)			
			TH-18	91 cm			FTHVS618..d-91 (Gore Reference: 6N16)	FTHVS518..d-91 (Gore Reference: 5N04)	FTHVS418..d-91 (Gore Reference: 4N02)	FTHVS318..d-91 (Gore Reference: 3N04)		
		22	TH-22	91 cm								FTHVS222..h-91 (Gore Reference: 2N03)
			22 (with pletgets)	TH-22	91 cm				FTHVS522..d-91xPL (Gore Reference: 3201*)	FTHVS422..d-91xPL (Gore Reference: 3200*)		

		Approx USP Size	7/0	6/0	5/0	4/0	3/0	2/0	0	
		Thread Diameter, mm	0,091-0,108	0,109-0,167	0,168-0,245	0,246-0,306	0,307-0,421	0,422-0,517	0,518-0,625	
		Gore Reference Thread Size:	CV-8	CV-7	CV-6	CV-5	CV-4	CV-3	CV-2	
		Needle Length, mm	Gore Reference Needle Code	Thread Length	Golnit Product Code (Gore Reference)					
Taper Point 1/2 Circle 	Double Armed 	26	TH-26	91 cm		FTHVS626..d-91	FTHVS526..d-91 (Gore Reference: 5N20)	FTHVS426..d-91 (Gore Reference: 4N06)	FTHVS326..d-91 (Gore Reference: 3N08)	FTHVS226..d-91 (Gore Reference: 2N02)
			TH-26	122 cm				FTHVS426..d-122 (Gore Reference: 4U06)	FTHVS326..d-122 (Gore Reference: 3U08)	FTHVS226..d-122 (Gore Reference: 2U02)
		30		91 cm					FTHVS330..d-91	
			TH-35; THX-36	91 cm				FTHVS435..d-91 (Gore Reference: 4N12)	FTHVS335..d-91 (Gore Reference: 3N06)	FTHVS236..hd-91 (Gore Reference: 3N08)
		35; 36	TH-35 THX-36	122 cm					FTHVS335..d-122 (Gore Reference: 3U06)	
Piercing Point 3/8 Circle 		9	PT-9	46 cm		FGDVS709..d-46 (Gore Reference: 7J04)				
			PT-9	61 cm	FGDVS809..d-61 (Gore Reference: 8K08)	FGDVS709..d-61 (Gore Reference: 7K04)	FGDVS609..d-61 (Gore Reference: 6K06)			
			PT-9	76 cm	FGDVS809..d-76 (Gore Reference: 8M08)	FGDVS709..d-76 (Gore Reference: 7M04)	FGDVS609..d-76 (Gore Reference: 6M06)			

		Approx USP Size	7/0	6/0	5/0	4/0	3/0	2/0	0			
		Thread Diameter, mm	0,091-0,108	0,109-0,167	0,168-0,245	0,246-0,306	0,307-0,421	0,422-0,517	0,518-0,625			
		Gore Reference Thread Size:	CV-8	CV-7	CV-6	CV-5	CV-4	CV-3	CV-2			
		Needle Length, mm	Gore Reference Needle Code	Thread Length	Golnit Product Code (Gore Reference)							
Piercing Point 3/8 Circle 	Double Armed 	13	PT-13	61 cm	FGDVS713..d-61 (Gore Reference: 7K08)	FGDVS613..d-61 (Gore Reference: 6K08)	FGDVS513..d-61 (Gore Reference: 5K08)					
			PT-13	76 cm	FGDVS813..d-76 (Gore Reference: 8M06)	FGDVS713..d-76 (Gore Reference: 7M08)	FGDVS613..d-76 (Gore Reference: 6M08)					
			PT-13	91 cm				FGDVS513..d-91 (Gore Reference: 5N08)				
		17	PT-17	91 cm				FGDVS517..d-91 (Gore Reference: 5N10)	FGDVS417..d-91 (Gore Reference: 4N10)	FGDVS317..d-91 (Gore Reference: 3N14)		
			Piercing Point 1/2 Circle 	13	PH-13	76 cm		FGHVS613..d-76 (Gore Reference: 6M14)	FGHVS513..d-76 (Gore Reference: 5M16)			
					PH-17	91 cm				FGHVS517..d-91 (Gore Reference: 5N14)	FGHVS417..d-91 (Gore Reference: 4N08)	FGHVS317..d-91 (Gore Reference: 3N12)
24, 25	PH-24	91 cm						FGHVS424..d-91 (Gore Reference: 4N16)	FGHVS324..d-91 (Gore Reference: 3N10)	FGHVS224..d-91 (Gore Reference: 2N10)		
	PH-24	122 cm						FGHVS424..d-122 (Gore Reference: 4U16)	FGHVS324..d-122 (Gore Reference: 3U10)	FGHVS224..d-122 (Gore Reference: 2U10)		
30	PH-30	107 cm							FGHVS330..d-107 (Gore Reference: 3P24)			
35, 36	PH-35	122 cm						FGHVS335..d-122 (Gore Reference: 3U20)				

Commonly Requested GOLNIT® PTFE Sutures

Proximal Aorto-bifemoral			
Approx USP Size (Gore Reference Thread Size):	Thread Length	Product Code (Gore Reference) with Taper Point Needles	Product Code (Gore Reference) with Tapercut Needles
2/0 (CV-3)	91 cm	FTHVS326..d-91 (Gore Reference: 3N08)	FGHVS324..d-91 (Gore Reference: 3N10)
3/0 (CV-4)	91 cm	FTHVS426..d-91 (Gore Reference: 4N06)	FGHVS424..d-91 (Gore Reference: 4N16)
	91 cm	FTHVS418..d-91 (Gore Reference: 4N02)	
Distal Aorto-bifemoral			
4/0 (CV-5)	76 cm	FTHVS513..d-76 (Gore Reference: 5M12)	
	76 cm	FTDVS513..d-76 (Gore Reference: 5M02)	
	91 cm	FTHVS518..d-91 (Gore Reference: 5N04)	
5/0 (CV-6)	76 cm	FTDVS612..d-76 (Gore Reference: 6M10)	FGDVS613..d-76 (Gore Reference: 6M08)
Chordae Tendinae			
3/0 (CV-4)	91 cm	FTHVS422..d-91xPL* (Gore Reference: 3200*)	
	91 cm	FTHVS418..d-91 (Gore Reference: 4N02)	
	122 cm	FTHVS418..-122 (Gore Reference: 4U02)	
4/0 (CV-5)	91 cm	FTHVS522..d-91xPL* (Gore Reference: 3201*)	
	91 cm	FTDVS417..d-91 (Gore Reference: 4N04)	
Proximal Femoro-popliteal			
4/0 (CV-5)	91 cm	FTDVS513..d-91 (Gore Reference: 5N02)	FGDVS513..d-91 (Gore Reference: 5N08)
5/0 (CV-6)	76 cm	FTDVS612..d-76 (Gore Reference: 6M10)	FGDVS613..d-76 (Gore Reference: 6M08)
Distal Femoro-popliteal			
5/0 (CV-6)	76 cm	FTDVS609..d-76 (Gore Reference: 6M02)	FGDVS609..d-76 (Gore Reference: 6M06)
6/0 (CV-7)	76 cm	FTDVS709..d-76 (Gore Reference: 7M02)	FGDVS709..d-76 (Gore Reference: 7M04)

* specification with ePTFE pledgets

Carotid Endarterectomy			
Approx USP Size (Gore Reference Thread Size):	Thread Length	Product Code (Gore Reference) with Taper Point Needles	Product Code (Gore Reference) with Tapercut Needles
5/0 (CV-6)	76 cm	FTDVS609..d-76 (Gore Reference: 6M02)	
	76 cm	FTDVS613..d-76 (Gore Reference: 6M04)	
Axillo-Femoral			
4/0 (CV-5)	91 cm	FTDVS513..d-91 (Gore Reference: 5N02)	FGDVS513..d-91 (Gore Reference: 5N08)
5/0 (CV-6)	76 cm	FTDVS609..d-76 (Gore Reference: 6M02)	
	76 cm	FTDVS612..d-76 (Gore Reference: 6M10)	
	76 cm	FTDVS613..d-76 (Gore Reference: 6M04)	FGDVS613..d-76 (Gore Reference: 6M08)
AV Access Graft			
5/0 (CV-6)	61 cm	FTDVS609..d-61 (Gore Reference: 6K02)	
	76 cm	FTDVS609..d-76 (Gore Reference: 6M02)	
	61 cm	FTDVS612..d-61 (Gore Reference: 6K10)	
	76 cm	FTDVS612..d-76 (Gore Reference: 6M10)	
	61 cm	FTDVS613..d-61 (Gore Reference: 6K04)	
	76 cm	FTDVS613..d-76 (Gore Reference: 6M04)	
6/0 (CV-7)	61 cm	FTDVS709..d-61 (Gore Reference: 7K02)	
	76 cm	FTDVS709..d-76 (Gore Reference: 7M02)	
AV Fistulas			
6/0 (CV-7)	46 cm	FTDVS709..d-46 (Gore Reference: 7J02)	
7/0 (CV-8)	46 cm	FTDVS809..d-46 (Gore Reference: 8J02)	

Orders:

+380 067 239 9197 (Ukraine), Viber, WhatsApp

+1 718 219 0731 (USA), Viber, WhatsApp

Email: info@golnit.com

Skype: GOLNIT LTD



Golnit Ltd

Kharkivske Shosse, 48

Kyiv, Ukraine 02160

ptfesutures.com

GOLNIT® is a trademark of Golnit Ltd.

GORE®, GORE-TEX® are trademarks of W. L. Gore & Associates.